






Laser Specifications:

- 4 Axis movement:
 - X axis - Horizontal movement:
20mm (0.8") +/-10mm (from ISO center)
 - Y axis - Vertical movement:
20mm (0.8") +/-10mm (from ISO center)
 - Roll axis - Horizontal rotation:
10° +/-5°
 - Pitch axis - Vertical rotation:
6° +/-3°
- Movement and position sensing detection for opposed lasers
- Crosshair beam rotation: 15° +/-7.5°
- Focus Adjustment: from 0.5m-6m (2'-20')
- Crosshair beam fan-out: 60° in both axis
- All displacements are precision motor driven.
- Laser Power Adjustment: 0.1mW-5mW
- Maximum line width @ 5m (16'): <0.5mm (0.02")
- Remote can control up to 7 lasers per room
- Operating Temperature: 5°-28°C (40°-85°F)
- Power requirements: 100-240 VAC, 0.25A, 50-60Hz
- Dimensions:
 - Length 210mm (8.3")
 - Width 116mm (4.5")
 - Height 94mm (3.7")
- Weight: 2.9 kg (6.4 lbs)



PATIENT POSITIONING LASER



Canada: 

Cemar Electro Inc.
1370 55th Ave.
Lachine, QC
Canada H8T 3J8

Tel: 514-631-5807
Fax: 514-631-7505
Toll Free: 1-800-298-5273
cemar@cemarelectro.com

United States: 

Cemar Electro Inc.
100 Walnut St.
Champlain, NY
USA 12919

Tel: 518-298-3065
Fax: 518-298-2927
Toll Free: 1-800-298-5273
cemarusa@cemarelectro.com

Europe: 

Scandinavia
<http://www.scanflex.se/>

Switzerland
<http://www.sirad.ch>

Coming Soon:
Spain UK
Portugal Russia
Ukraine Middle East


WWW.CEMARELECTRO.COM





- Alignment verification sensor
- Auto ambient light adjustment
- 4 Axis movement
- Precision motor drives
- Micro and macro adjustment
- Available in **Green**, **Red**, and **Blue**


Optional Mounting Bracket:




Remote control touch-screen adjustments:



Main Menu



Adjustment Selection



Continuous Adjustment

- Touch screen remote
- Full range, non-directional
- User friendly icons
- Robust ergonomic design
- RF wireless
- Full function

