

MacroMedics®

Vacuum Cushions

MacroMedics® easy-to-use vacuum cushions are manufactured of first class, high quality, durable and easy to clean materials. The vacuum cushion is a well known and widely accepted aid for patient immobilisation in radiotherapy.

MacroMedics vacuum cushions are filled with small polystyrene pellets to create a firm but comfortable cradle to accurately support and position the patient. They naturally adapt to the patients anatomical contours. A wide range of dimensions and filling volumes are available.

The vacuum cushions feature:

- Easy-to-use, reusable vacuum cushions
- Filled with small polystyrene pellets to create a firm but comfortable cradle around the patient
- The vacuum cushions are radiotranslucent and easy to clean
- Durable, high quality blue Polyamide cushions
- Most common types quick-release valves are available

155110	VHS3040	Vacuum Cushion, head support, 30x40cm
155140	VHS2550	Vacuum Cushion, head support, 25x50cm
155150	VHS5070	Vacuum Cushion, head & shoulder support, T-shape, 50x70cm
155210	VBR5070	Vacuum Cushion, breast support, 50x70cm
155300	VPS6565	Vacuum Cushion, pelvic support, 65x65cm
155400	VBOS10070	Vacuum Cushion, body support, 100x70cm
155410	VBOS10080	Vacuum Cushion, body support, 100x80cm
155430	VBOS12080	Vacuum Cushion, body support, 120x80cm
155450	VBOS100150	Vacuum Cushion, body support, 100x150cm
155510	VBOS200100	Vacuum Cushion, body support, 200x100cm

Contact MacroMedics or your local distributor for detailed information or an on-site product demonstration.


macromedics

SOLUTIONS IN RADIOTHERAPY

MacroMedics BV
Kouwe Hoek 18
2741 PX Waddinxveen
The Netherlands

Tel: +31(0)182 389777
Fax: +31(0)182 389778
info@macromedics.com
www.macromedics.com